

InterSystems Official
[Mike Morrissey](#) · 20 mayo, 2020

May 20, 2020 - Alert: HS2020-05: HealthShare Alerts

Dear HealthShare Customer:

This post is part of the HealthShare HS2020-05 Alert communications process. The same information is also distributed:

- By Email
- On the [Product Alerts & Advisories page](#)
- On the [WRC Distribution Page InterSystems Documents](#)

There are 3 alerts in the HealthShare HS2020-05 Alert communication, including the previously posted "[Alert: Possible Data Integrity Issues with Online Backup of Large Databases](#)". The Alert Summary is in the table below, and the detail is contained in the attached document: [HS2020-05-Communication](#)

Alert	Product & Versions Affected	Risk Category & Score
HS2020-05-01: On an Android device that is in sleep mode, the standalone v2 Clinical Viewer using Chrome does not time out	The affected products & versions are: <ul style="list-style-type: none">• Information Exchange 2018.1• Unified Care Record 2019.1, 2019.1.2, 2019.2	3-Medium Risk (Security)
HS2020-05-02: FHIR Request from ODS Not Logged in ATNA if Custom Pairs are Used on the Patient Streamlet	The affected products & versions are: <ul style="list-style-type: none">• Information Exchange 2018.1• Unified Care Record 2019.1, 2019.2, 2020.1	3-Medium Risk (Privacy)
HS2020-05-03: Possible Data Integrity Issues with Online Backup of Large Databases	The affected products & versions are: <ul style="list-style-type: none">• All versions of all InterSystems HealthShare products	5-Very High Risk (Operational)

If you have any questions regarding this advisory, please contact the Worldwide Response Center (WRC) at support@InterSystems.com or +1.617.621.0700., and reference " HealthShare Alert HS2020-05 " .

[#Alertas](#) [#HealthShare](#) [#InterSystems Official](#)

URL de fuente: <https://es.community.intersystems.com/node/477176>