
Pregunta

[Michel Liberado](#) · 1 mar, 2020

Unable to retrieve response object allocated in the last method of a nested SendRequestSync sequence

Hi,


I have a production with a BS calling a BP calling another BP and so on and finally a BO. All calls are made using SendRequestSync.

The initial request is passed by value. The response is passed by ref and I allocate memory for it in the BO, hoping that it will go back to the initial caller in the BS, but it does not.

It looks like the following:

```
// In BS, do not allocate response, call BP1
Set tSC = SendRequestSync(request, .response)
// In BP1, call BP2. Now I don't know if .response is a reference to a reference
Set tSC = SendRequestSync(request, .response)
// In BP2, call BP3. Now I don't know if .response is a reference to a reference to a reference
Set tSC = SendRequestSync(request, .response)
// Do some stuff with reqa and reqb an then the last one supposed to allocate response
Set tSC = SendRequestSync(reqc, .response)
```

I have tried many combinations of reference of values, does not work. Why is the response NULLED and I can't retrieve it in the initial caller ?



#Operación empresarial #InterSystems IRIS for Health

URL de fuente: <https://es.community.intersystems.com/node/474456>